

Critical and Volumetric Flow Controllers from FIA

Control and monitoring of the air-flow is pivotal when it comes to performance testing of inhalation devices. FIA offers three models ranging from our standard trigger box to our most advanced critical flow controller.

Typical usages are for delivered dose uniformity and aerodynamic particle size distribution testing. The products have been designed from both the user and regulatory perspective, why the equipment is compact and user-friendly but still complies with the European Pharmacopoeia and the US Pharmacopeia. Acknowledging a lab set up may become crowded with lots of tubing and electrical cords, FIA's range of equipment minimises the number of units the analyst needs to handle on the bench.



Trigger Model I



Optional output and inputs configurable according to customer needs.

Trigger Box Model I – Standard

This model gives the analyst the basic functionality in a very compact and easily managed format.

Key features are:

- Generous flow path which minimizes loss of flow and helps attaining the critical flow condition (P3/P2 ratio)
- Flow control valve
- Timer controlled valve
- Resettable counter for actuations
- Optional
 - Input to actuate from external equipment or from optional handheld or foot switch
 - Configurable electrical, logical and pneumatic output to actuate other equipment (such as a camera) at a set time after start.





www.fia.se



www.fia.se



Trigger Box Model II – Critical Flow Controller

This model builds on the functions of Model I with the addition of differential pressure transducers for measuring and presenting pressure at P1, P2 and P3.

Key features are:

- Generous flow path which minimizes loss of flow and helps attaining the critical flow condition
- Flow control valve
- Timer controlled valve
- Measurement of P1 and critical flow (P3/P2 ratio)
- Resettable counter for actuations
- 3.5" TFT touch screen for control and presentation of measured values, suitable to operate with or without gloves
- Port for external printer (list of dose #, P1 and ratio)
- Input to actuate from optional handheld or foot switch

For contact and more information

Kjell Fransson AB FIA

Managing Director kjell.fransson@fia.se Mobile: +46 701 45 54 83





AB FIA Vinkelhaken 1 D SE-247 32 Södra Sandby, Sweden +46 46 13 90 00 i@fia.se • www.fia.se



Trigger Box Model III – Critical Flow Controller with Flow Measurement

Trigger Box Model III has the same functions as Model II, but is also equipped with an internal flow meter. With this equipment the user has full control of the testing, including volumetric flow measurement. The latter is enabled by an integrated laminar flow element downstream the test set-up (e.g. dose collector or impactor). By using a conversion factor set at the factory and simultaneous pressure measurement, the air flow into the device can be calculated. The principle has proven to be robust and has been used by FIA's customers for over 25 years but has now been integrated into a compact format for routine use.

Key features are:

- Generous flow path which minimizes loss of flow and helps attaining the critical flow condition
- Flow control valve
- Timer controlled valve
- Measurement of P1 and critical flow (P3/P2 ratio)
- · Measurement of volumetric air flow into the device
- Resettable counter for actuations
- 7" TFT touch screen for control and presentation of measured values, suitable to operate with or without gloves
- Ports for external printer (list of dose #, P1, ratio and Q) and USB for export to usb-stick in .csv format
- Input to actuate from optional handheld or foot switch